

Medit Link & Apps

是一個整合診所與牙體技術所工作流程的開放性系統平台



分析圖表



工作流程管理



雲端系統



病患溝通

以病患為中心的工作流程
與合作夥伴共同創造價值
人性化的數據管理系統

透過 Medit Link 體驗數位牙科的魅力

掃描軟體具備免費的智慧功能與應用程式
協助醫師打造高品質精準的治療方案
提供病患舒適良好的看診體驗




臨時假牙



建立模型



內冠鬆緊度



影像整合準確度



微笑設計



DCM 轉檔



咬合分析



矯正模擬



校正工具



www.medit.com

欲了解更多 Medit i700 相關資訊及操作方式
請掃描下方 QR code.



User Group



Medit Academy



Help Center



Contact Medit

智能 精準 快速 MEDIT i700



DIO IMPLANT

台北總公司 TEL: (02) 8752-3319
台南辦事處 TEL: (06) 302-2521

台中辦事處 TEL: (04) 2305-0800
高雄辦事處 TEL: (07) 550-1859

DIO IMPLANT

Medit i700 提升牙科治療與診所效率的一大關鍵！

- 流暢、人性化的外觀設計
- 智能軟體 Medit Link 整合所有工作流程

取代傳統繁複的印模過程，智能、精準、快速的 i700，為醫師與病患提供優質的掃描體驗



reddot winner 2021



245g
輕量外型



Up to 2X
極致快速



UV-C LED
紫外線殺菌



180°
旋轉式掃描頭



遠端
操控模式



拆卸式
電源線

規格 Specifications

Category		Medit i700
Scanning technology	Scanning frame	Up to 70 FPS
	Imaging technology	3D-in-motion video technology 3D full color streaming capture
	Light source	LED
Accuracy	Anti-fogging technology	Adaptive anti-fogging
	Full-arch	10.9µm ± 0.98
Handpiece	Dimensions	248 x 44 x 47.4 mm
	Weight	245 g
Tip	Tip size	22.2 x 15.9 mm
	Mirror angle	45-degree angle <small>Easier to scan distal molar area</small>
	Scan area	15 x 13 mm
	Autoclavable	Up to 100 times Autoclave 121°C 30 mins Autoclave 134°C 4 mins
	Reversible tip	Yes
Special mode	Remote control mode	Yes
	UV-C disinfection	Yes
Cable	Length	2.0 m / Handpiece Detachable <small>2.5m (Sold Separately)</small>
	Connectivity	USB 3.1 Gen1 (C Power Delivery)

請掃描右方 QR code
以了解系統需求的詳細資訊 >>

